FA AVAILABILITY & PRICING

Product Number	Product Name	Description	Packaging	Pricing	Competitive Products
34471	HYDRANAL NEXTGEN Coulomat A-FA (anolyte)	Anolyte for coulometric Karl Fischer titration of ketones and Li-ion battery electrolytes, preferred for cells with diaphragm, acetonitrile-based, free of alcohols (free of CMR substances)	500 ml	+23% over Coulomat AK	No competitive products!
34470	HYDRANAL NEXTGEN Coulomat C-FA (catholyte)	Catholyte for coulometric Karl Fischer titration of ketones and Li-ion battery electrolytes, acetonitrile-based, free of alcohols (free of CMR substances)	10 × 5 ml ampoules	+29% over Coulomat CG-K	

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
TITRATION PROBLEMS WITH CURRENTLY AVAILABLE KF REAGENTS

- Alcohols undergo side reactions with the additives used in many LiB electrolytes.
- Difficult or even impossible to measure water content in most prominent LiB electrolytes that contain additives such as:
 - VC
 - FEC
 - Borates such as LiBOB and LiDFOB

Sample	Coulomat A-FA + C-FA	Coulomat AK + CG-K	Reagent M
LiB electrolyte with VC/FEC	+	o	o
LiB electrolyte with LiBOB	+	-	-
LiB electrolyte with LiDFOB	+	-	-
Pure VC	+	-	-
Pure FEC	+	o	o
Pure LiBOB	+	-	-
Pure LiDFOB	+	-	-
Ketones	+	+	-

+ accurate titration, suppressed side effects
 o limited titration (possible only for small amount of sample)
 - no accurate titration possible

VC = vinylene carbonate
 FEC = fluoroethylene carbonate
 LiBOB = lithium bis(oxalate)borate
 LiDFOB = lithium difluoro(oxalato)borate

Honeywell Confidential - ©2022 by Honeywell International Inc. All rights reserved.

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SIDE REACTIONS: UNSATURATED COMPOUNDS

Sample	As raw material	As additive	Comment
Unsaturated compounds (e.g. VC, vinyl acetate)	No or fading end point Accurate water determination is impossible.	Increasing drift. Decreasing accuracy of titration.	Iodine-consuming side reaction occurs.
FEC	Causes fewer problems compared to VC when sample size is small Otherwise, no or fading end point.	Causes fewer problems if the absolute water content in the sample is low.	Can be titrated with low water content in small portions; however, small sample size increases weighing error.


Fig. 2. Titration of pure VC.

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


Fig. 2. Titration of pure VC.

SIDE REACTIONS: BORATES

Sample	As raw material	As additive	Comment
Borates (e.g. LiBOB, LiFOB)	Accurate water determination is impossible.	Long titration with highly overdetermined results. Note: end point can be reached, but giving incorrect result	Water-releasing side reaction occurs.

Honeywell Confidential - ©2022 by Honeywell International Inc. All rights reserved.

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1. Titration of 1 ml 5 wt% LiBOB in acetonitrile

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

SIDE REACTIONS: KETONES

Sample	As raw material	As additive	Comment
Ketones (e.g. acetone)	No or fading end point. Accurate water determination is impossible.	Long titration with highly overestimated results. No or fading end point. Increasing drift.	Water-releasing side reaction occurs. Acetone and diketones are decomposition products in LIB electrolytes containing propylene carbonate (PC).

Honeywell Confidential - ©2022 by Honeywell International Inc. All rights reserved.

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